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
Fig. 5. Topographic plot of GC-IMS spectra. (A) Pomace olive oil, (B) olive oil and (C) extra virgin olive oil.


Fig. 3. PCA results for mean-centered data for the three types of oil as obtained with UV-IMS. Samples: Pomace olive oil (red triangles); olive oil (green asterisks); extra virgin olive oil (blue squares). (For interpretation of the references to color in this figure legend, the reader is referred to the web version of the article.)

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Reference: Direct classification of olive oils by using two types of ion mobility spectrometers. Anal. Chim.Acta (2011) .DOI: 10.1016/j.aca.2011.03.007.

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Reference: Ion mobility spectrometry coupled to gas chromatography: A rapid tool to assess eggs freshness. **Food Chemistry**, 271 (2019) 691-696.

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Reference: Target vs spectral fingerprint data analysis of Iberian ham samples for avoiding labelling fraud using headspace–gas chromatography–ion mobility spectrometry. **Food Chemistry**, 2. 46(2018) 65-73.

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Reference: Characteristic volatiles fingerprints and changes of volatile compounds in fresh and dried *Tricholoma matsutake* Singer by HS-GC-IMS and HS-SPME-GC-MS. *Journal of Chromatography B*, 1099. (2018) 46-55.

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Referenc: Volatile-Compound Fingerprinting by Headspace-Gas-Chromatography Ion-Mobility Spectrometry (HS-GC-IMS) as aBenchtop Alternative to 1H NMR Profiling for Assessment of the Authenticity of honey. *Analytical Chemistry* DOI: 10.1021/acs.analchem.7b03748 (2017)

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Reference: Content variations in compositions and volatile component in jujube fruits during the blacking process. **Food Sci Nutr.** 2019;1-9.

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- 1.A green triple-locked strategy based on volatile-compound imaging, chemometrics, and markers to discriminate winter honey and sapium honey using headspace gas chromatography-ion mobility spectrometry[J]. **Food Research International**, <https://doi.org/10.1016/j.foodres.2019.01.004>.
- 2.Development of a flavor fingerprint by HS-GC-IMS with PCA for volatile compounds of Tricholoma matsutake Singer. **Food Chemistry**, 290 (2019) 32-39.
- 3.Characterization of the aroma release and perception of white bread during oral processing by gas chromatography-ion mobility spectrometry and temporal dominance of sensations analysis. **Food Research International**, 123 (2019) 612-622.
- 4.robustness study of calibration models for olive oil classification: Targeted and non-targeted fingerprint approaches based on GC-IMS. **Food Chemistry**, 288 (2019) 315-324.
- 5.Direct authentication of three Chinese materia medica species of the Lili Bulbus family in terms of volatile components by headspace-gas chromatography-ion mobility spectrometry. **Analytical Methods**, 2019. DOI: 10.1039/c8ay02338g.
- 6.Content variations in compositions and volatile component in jujube fruits during the blacking process. **Food Sci Nutr.** 2019. 1-9.
- 7.Identification of terpenes and essential oils by means of static headspace gas chromatography-ion mobility spectrometry. **Anal Bioanal Chem**, 409 (2017) 6595-6603.
- 8.Detection of Adulteration in Canola Oil by Using GC-IMS and Chemometric Analysis. **International Journal of Analytical Chemistry**, 2018, Article ID 3160265.
- 9.Characteristic volatiles fingerprints and changes of volatile compounds in fresh and dried Tricholoma matsutake Singer by HS-GC-IMS and HS-SPME-GC-MS. **Journal of Chromatography B**, 1099 (2018) 46-55.
- 10.Ion mobility spectrometry coupled to gas chromatography: A rapid tool to assess eggs freshness. **Food Chemistry**, 271 (2019) 691-696.
- 11.Volatile-Compound Fingerprinting by Headspace-Gas Chromatography Ion-Mobility Spectrometry (HS-GC-IMS) as a Benchtop Alternative to 1H NMR Profiling for Assessment of the Authenticity of Honey. **Analytical Chemistry**, DOI: 10.1021/acs.analchem.7b03748.
- 12.Coupling Gas Chromatography with Ion Mobility Spectrometry. **LC*GC Europe**, 2016. 294-303.
- 13.Target identification of volatile metabolites to allow the differentiation of lactic acid bacteria by gas chromatography-ion mobility spectrometry. **Food Chemistry**, 220 (2017) 362-370.

14. Ion mobility spectrometry of volatile compounds from Iberian pig fat for fast feeding regime authentication. **Talanta**, **76** (2008) 591–596.
15. A review of recent, unconventional applications of ion mobility spectrometry (IMS). **Analytica Chimica Acta**, **703** (2011) 114–123.
16. Using ion mobility spectrometry for screening the autoxidation of peanuts. **Food Control**, **64** (2016) 17–21.
17. Ion mobility spectrometry fingerprints: A rapid detection technology for adulteration of sesame oil. **Food Chemistry**, **192** (2016) 60–66.
18. Molecular Analysis for Determining Differences and Similarities of Perfume Dupes. **GLOBAL INGREDIENTS & FORMULATIONS GUIDE 2016**.
19. Ion Mobility Spectrometry versus Classical Physico-chemical Analysis for Assessing the Shelf Life of Extra Virgin Olive Oil According to Container Type and Storage Conditions. **J. Agric. Food Chem**, **63** (2015) 2179–2188.
20. Gas chromatography and associated techniques in the analysis of laundry malodour. **H&PC Today - Household and Personal Care Today**, Vol. 9 nr.
21. Evaluation and application of static headspace–multicapillary column–gas chromatography–ion mobility spectrometry for complex sample analysis. **Journal of Chromatography A**, **1338** (2014) 136–148.
22. Direct determination of 2,4,6-trichloroanisole in wines by single-drop ionic liquid microextraction coupled with multicapillary column separation and ion mobility spectrometry detection. **Journal of Chromatography A**, DOI:10.1016/j.chroma.2011.06.032.
23. Volatile Metabolites of Goat Cheeses Determined by Ion Mobility Spectrometry. Potential Applications in Quality Control. **Food Anal Methods**, DOI: 10.1007/s12161-014-0050-1.
24. Multi-capillary column-ion mobility spectrometry: a potential screening system to differentiate virgin olive oils. **Anal Bioanal Chem**, **402** (2012) 489–498.
25. Ion Mobility Spectrometry a Versatile Analytical Tool for Metabolomics Applications in Food Science. **Current Metabolomics**, **2** (2014) 264–271.
26. Direct coupling of a gas–liquid separator to an ion mobility spectrometer for the classification of different white wines using chemometrics tools. **Talanta**, **84** (2011) 471–479.
27. Determination of volatile compounds by GC–IMS to assign the quality of virgin olive oil. **Food Chemistry**, **187** (2015) 572–579.
28. Resolution-optimized headspace gas chromatography-ion mobility spectrometry (HS-GC-IMS) for non-targeted olive oil profiling. **Anal Bioanal Chem**, DOI: 10.1007/s00216-017-0338-2.

29. Direct classification of olive oils by using two types of ion mobility spectrometers. **Analytica Chimica Acta**, DOI:10.1016/j.aca.2011.03.007.
30. Target vs spectral fingerprint data analysis of Iberian ham samples for avoiding labelling fraud using headspace – gas chromatography–ion mobility spectrometry. **Food Chemistry**, **246** (2018) 65–73.
31. Direct differentiation of herbal medicine for volatile components by a multicapillary column with ion mobility spectrometry method. **J. Sep. Sci**, **38** (2015) 3205–3208.
32. Use of a non-destructive sampling method for characterization of Iberian cured ham breed and feeding regime using GC-IMS. **Meat Science**, <https://doi.org/10.1016/j.meatsci.2019.02.018>.
33. Biochemical changes induced by dominant bacteria in chill-stored silver carp (*Hypophthalmichthys molitrix*) and GC-IMS identification of volatile organic compounds. **Food Microbiology**, **84** (2019) 103248.
34. HS-GC-IMS and chemometric data treatment for food authenticity assessment: Olive oil mapping and classification through two different devices as an example. **Food Control**, DOI: <https://doi.org/10.1016/j.foodcont.2018.11.001>.
35. Quality assessment of olive oils based on temperature-ramped HS-GC-IMS and sensory evaluation: Comparison of different processing approaches by LDA, kNN, and SVM. **Food Chemistry**, DOI: <https://doi.org/10.1016/j.foodchem.2018.11.095>.

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